

MobyNeb

Dual-purpose
aerosol therapy
device



Mobyneb features the new Perfecta spacer mask for adults and children, as well as Rinowash, which can be used

to treat localized disorders of the upper and lower airways. It is the first device to incorporate the innovative concept

of DISTRICT INHALATION THERAPY, which facilitates drug deposition in the target organ.




Mobyneb is comprised of the following:

- Power unit with accessory compartment
- Perfecta
- Rinowash nebuliser
- Connecting tube



Innovative spacer mask for adults and children, complete with Mefar2000 nebuliser. Specific for treating disorders of the lower airways. Allows high deposition of the atomised drug.



Used together with the **Rinowash** nasal nebuliser, Mobyneb allows physiological and/or medical flushing of the nasal cavities. Rinowash produces a jet of atomised solution, which can be directed into the nasal cavities to enhance hydration and fluidification of the mucus and catarrh, and hence their removal. It comes with a single bell suitable for all ages.



N.B.: Perfecta and Rinowash are also sold separately and can be connected to any pneumatic aerosol therapy device.



As from 1st June 2009 Markos Mefar has become Air Liquide Medical Systems S.p.A.

AIR LIQUIDE MEDICAL SYSTEMS S.p.A.
Via dei Prati 62 - 25073 Bovezzo (Brescia)
Tel. +39 030 20159.11
Fax +39 030 2098329

www.airliquidemedicalsystems.it
CE₀₀₅₁ www.airliquidemedicalsystems.com

Air Liquide's Healthcare business supplies medical gases & equipment to hospitals, medical services to patients at home and is active in the fight against nosocomial infections. Its 7,800 employees serve 5,000 hospitals and 300,000 patients, on five continents. Air Liquide is the world leader in gases for industry, health and the environment.